


Single-cell methodologies and the ESCG facility

SciLifeLab

Karolina Wallenborg
January 25, 2021




Eukaryotic Single-cell Genomics (ESCG)

- Established in 2015 with funding from SciLifeLab
- 2021
 - 6 personnel




Rickard Sandberg Sten Linnarsson



Anastasios, Karolina, Matilda, Henrik, Marcela, Samaneh & Helena




Fluidigm C1 Auto-prep single cell system




Aldridge & Teichmann, Nature Communication, 2020

Single cell RNA-sequencing experiment




- User Consultations
 - Experimental design
- Quality control of single-cell preparations prior to the experiment
 - Sample submission guidelines for cells/nuclei for different applications
- Initial bioinformatics

Single-cell isolation and library preparation



Aldridge & Teichmann, Nature Communications, 2020



RNA molecule capture and reverse transcription

Amplification

Library preparation

3' -tag

5' - tag



Full-length

Adapted from
Kolodziejczyk A et al, Molecular Cell, 2015

Full-length vs 3' or 5'-end transcript sequencing


- Full-length
 - Whole transcript information
 - Gene expression quantification
 - Isoform, SNPs, and mutations
 - Higher sensitivity
- Tag-based methods
 - Estimate of transcript abundance
 - Early multiplexing (lower cost)
 - Combined with molecular counting
 - Retain DNA strand information
 - Higher throughput

Full-length transcriptome coverage + 5' UMI RNA counting with Smart-seq3



- Highly sensitive (80% of molecules detected by smRNA-FISH per cell)


Single-cell ATAC-seq (Assay for Transposase Accessible Chromatin)



(Buenrostro et al, Nature Methods, 2013)


Chromatin and DNA-binding proteins regulate gene expression

- Subpopulations with different chromatin accessibility profiles
→ Increased understanding of gene regulatory networks upstream on gene expression




Adapted from Chen et al, Genome Biology, 2019

Single-cell Immune Profiling, V(D)J



Every
receptor is
unique


Example from 10XGenomics protocol




cDNA amplification

Target enrichment with specific primers
for TCR or Ig constant regions

Chromium Single Cell V(D)J Enriched Library




Chromium Single Cell 5' Gene Expression Library




CITE-seq

simultaneously measure gene expression and cell surface protein abundance



Adapted from Stoeckius et al., 2017




Stoeckius et al., 2017

Antibody panels can
be purchased from
Biolegend

→ Integrated protein and RNA measurements enhance cellular phenotyping

Cell hashing – multiplexing with ubiquitous surface protein expression




- Cell hashing, Stoeckius, Zheng et al., 2017
- Sample multiplexing
- Identify doublets, superload
- Batch effects
- Save on costs

Sample multiplexing

- Demuxlet (genetic variation), Kang et al, Nature Biotech, 2018
- MULTI-seq, McGinnis et al Nature Methods, 2019
 - Lipid-tagged indices
- Cell hashing, Stoeckius et al Genome Biology, 2018
- Nuclei multiplexing with barcoded antibodies
- CellTag indexing, Guo et al, Genome Biology, 2019
- Multiplexing by transient barcoding, Shin et al Science Advances, 2019

Single-cell multimodal omics



scRNA-seq + scATAC-seq
(SNARE-seq)

scCUT&Tag

Zhu et al, Comment in Nature Methods, 2020


Services at ESCG

- Single-cell RNAseq
- Single-nucleus RNAseq
- Single-cell ATACseq
- Single-cell CUT&Tag
- Single-cell Immune profiling
- CITE-seq and cell/nuclei hashing
- Single-cell CRISPR screens
- Sequencing is done with NGI and Clinical Genomics in Stockholm

NEW: 10XGenomics Chromium platform in BSL3 for infections biology samples.


Services and applications at ESCG

Gene expression




- 10XGenomics
 - 3'GE
 - 5'GE
- Smart-seq2
- Smart-seq3

Gene expression + Immune profiling




- 10X + VDJ
 - 5'GE
 - CITE-seq
 - Cell hashing

Gene expression + Cell surface protein expression



- 10X + CITE-seq
 - 3' GE
 - 5' GE
- Cell hashing
 - 3' GE
 - 5' GE


Gene expression + CRISPR screen



- 10X + CRISPR
 - 10X 3' GE

High-throughput
Genome engineering
facility

Chromatin accessibility



- 10XGenomics
 - scATAC-seq


ESCG should enable scale in single-cell genomics hard to achieve in individual labs



TECAN EVO, liquid handling robot




ICELL8cx, TAKARA



9600 wells



5184 wells



NanoDrop dispenser, GC Biotech



Chromium Controller, 10XGenomics



8 x 10,000 cells

Outlook for single-cell genomics

- Multi-omics -- combine scRNA-sequencing with:
 - Surface protein expression (CITE-seq, Cell hashing)
 - CRISPR perturbations
 - scRNAseq + scATAC-seq
 - Histone modifications (scChIC-seq, CUT&Tag)
 - Low-coverage DNA (DLP+, Vancouver)
 - Methylated DNA (Joe Ecker lab, SALK)
- Spatial techniques
- Improve quality and reduce costs

Sample preparation

- Limit time of cell isolation
 - Be gentle
 - Viability >80%
 - Careful cell handling
 - No debris in the suspension
 - No cell aggregations
-
- Nuclei
 - Cryo-preserved cells
 - Fixed cells (methanol)



Contact us for tips & tricks!

Contacts



**Karolinska
Institutet**

SciLifeLab

escg@scilifelab.se

<https://www.scilifelab.se/facilities/eukaryotic-single-cell-genomics>

<https://ki.se/en/research/core-facilities-for-research>

→ iLAB login

Our publications can be found here:

[https://publications.scilifelab.se/label/Eukaryotic%20Single%20Cell%20Genomics%20\(ESCG\)](https://publications.scilifelab.se/label/Eukaryotic%20Single%20Cell%20Genomics%20(ESCG))